

Thin Films Nanolithography (XEL)



Sample:
Nanolithography on
Silicon Oxide

Image Conditions:
Scan Size: (10um)
Pixel Resolution
(256x256)
Writing Voltage (5V)
Writing Mode (Raster)

System Requirement:
Closed Loop XY
Scanner, XEL,
Tip/Sample Bias


The Benefits

The XE-series AFMs have a very user friendly design, which allows users to customize the system for their advanced applications. One of the many options available is XEL. XEL is a nanolithography and nano-manipulation package for the XE-series AFM. It is a very easy to use software option with load, scanner position, and electrical bias patterning capability.


Park Systems Inc
3040 Olcott St
Santa Clara, CA 95054
Tel: 408.986.1110
Fax: 408.986.1199
www.parkAFM.com
info@parkAFM.com

Scanning Probe Nanolithography is a very promising technology for nanofabrication. It uses the probe tip of scanning probe microscope to manipulate materials at nanometer scale. There are many ways to make patterns, such as mechanical scratching and electrical field oxidation. Park System provides a very powerful software, XEL, for Nanolithography. It supports a variety of lithography modes such

Figure 1. A Big Buck pattern made in a silicon substrate using bias mode XEL Nanolithography. The pattern is made on a silicon surface using the Raster scan mode.



Publications Using XE-AFM and Nanolithography



AA Milner, K Zhang, and Y Prior; Floating Tip Nanolithography; NanoLetters; 2008 Vol. 8, No. 7; 2017-2022

Equipment: Park Systems XE-120

Figure 2. (a) Hot floating tip writing WIS (for the Weizmann Institute of Science) on AZ4620 photoresist film. The writing speed was 50 nm/sec at tip-sample gap 3 (2 nm and the average laser power was $\sim 3 \times 10^5$ W/cm². (b) Direct mechanical scratching of the same polymer with an identical cold tip. Profiles c and e correspond to the red cross-section lines in panels a and b; images d and f are the magnified line intersections in panels a and b, respectively.

Abstract: We demonstrate noncontact, high quality surface modification of soft and hard materials with spatial resolution of ~ 20 nm. The nanowriting is based on the interaction between the surface and the tip of a standard atomic force microscope illuminated by a focused femtosecond laser beam and hovering (at ambient conditions) 1-4 nanometers above the surface without touching it. Field enhancement at the tip-sample gap or high tip temperature are identified as the causes of material ablation.